

POX-OEM Serial Communication Protocol Revision 1.1

Mediaid Inc,
#17517 Fabrica Way Suite H,
Cerritos, CA 90703
USA

1. Scope

The purpose of this document is to define the serial communication protocol between all Mediaid POX OEM boards and OEM equipment.

2. Electrical/Timing

A three wire RS-232 style interface is all that is needed to communicate with the POX-OEM; RxD, TxD and ground. The RxD carries serial data to the POX_OEM while the TxD carries serial data from the POX_OEM. These two signals are TTL level signals. A TTL high is an RS-232 "mark" and a TTL low is an RS-232 "space". The characters are transferred at a rate of 19200 baud.

3. Command/Response packet format

The character format is: one start bit, 8 data bits (lsbit first) one stop bit. Multiple characters make up a command/response packet. The N+2 characters in the command/response packet are: "command"(or "response"), "data1", "data2",..., "dataN", "checksum" N is determined by the command/response and may be zero.

The timing between the "command/response" and "checksum" must not exceed two hundred fifty milliseconds (250ms).

The command characters (information from the OEM product to the POX-OEM) are the ASCII characters from "!" to "?". The response characters (information from the POX_OEM to the OEM product) are the ASCII characters from "a" to "~". They are all defined in the following section.

All dataN used in the command packets are base 32 digits where the base 32 zero character is ASCII "@", the base 32 thirty-one character is ASCII "_". This makes all characters displayable and contiguous. If a number to be transmitted is greater than 31 then it is split into groups of 5 bits and sent most significant first. Example: the number 263 is split as an 8 and 7 $[(8*32) + 7 = 263]$, the 8 is sent first.

The checksum is the negative modulo 32 summation of all the previous characters in the packet (including the command/response) plus the offset of 0x41. That is if all characters are modulo 32 summed in a packet their sum would be zero. (Hint: To create a checksum add all transmitted characters, negate, AND that negated sum with 0x1f and add 0x41, makes all checksum characters "@" through "_", same as the dataN characters.)

A graphical representation is:

Error data can range from 0 to 1023. There may be multiple errors in a response.

4.2.5.1. Error codes

4.2.5.1.1. No error, 0

4.2.5.1.2. ROM checksum error, 1(0x001)

4.2.5.1.3. Low power supply, 2(0x002)

4.2.5.1.4. EEPROM error, 4(0x004)

4.2.5.1.5. No red LED, 8(0x008)

4.2.5.1.6. No Ir LED, 16(0x010)

4.2.5.1.7. Thin tissue, 32(0x020)

4.2.5.1.8. Thick tissue, 64(0x040)

4.2.5.1.9. Maximum perfusion, 128(0x080)

4.2.5.1.10. System failure error, 256(0x100)

4.2.5.1.11. No module attached, 512(0x200)

4.2.5.1.12. Analog output calibration failure, 1024(0x400)

4.2.6. ASCII f Parametric (incomplete, x = data, chk = checksum)

```

f x x x x x x chk
| | | | | | | | _____ modulo 32 checksum
| | | | | | | | _____ red pulse height gain
| | | | | | | | _____ red perfusion gain
| | | | | | | | _____ red LED drive
| | | | | | | | _____ IR pulse height gain
| | | | | | | | _____ IR perfusion gain
| | | | | | | | _____ IR LED drive
| | | | | | | | _____ "parametric" character

```

4.2.7. ASCII g Version (incomplete, x = data, chk = checksum)

```

g x x x chk
| | | | |_____ modulo 32 checksum
| | | | |_____ software version (3rd byte - lsb)
| | | | |_____ software version (2nd byte)
| | | | |_____ software version (1st byte - msb)
| | | | |_____ "s/w version" character

```

The software version number can range from 0 to 32,767.

4.2.8. ASCII i Model number (incomplete, x = data, chk = checksum)

```

i x x x chk
| | | | |_____ modulo 32 checksum
| | | | |_____ model number (lsb)
| | | | |_____ model number
| | | | |_____ model number (msb)
| | | | |_____ "model number" character

```

The model number can range from 0 to 32,767.

4.2.9. ASCII j NAK (complete)

j@V NAK because of bad command
jAU NAK because of check sum
jBT NAK because of internal error
jCS NAK because of time out
jDR NAK because of bad parameter

These NAK commands (Not Acknowledge) are transmitted upon there being an error in the last OEM data packet, or, an internal POX error.

4.2.10. ASCII k ACK (complete)

kU

This ACK command (Acknowledge) is transmitted upon acceptance of the last transmitted command.

4.2.11. ASCII l Sensor type (incomplete, x = data, chk = checksum)

```

l x chk
| | | | |_____ modulo 32 checksum
| | | | |_____ sensor type
| | | | |_____ "sensor type" character

```

The sensor type can range from 0 to 31.

1 = M120 finger sensor module
2 = M120 cable adaptor module
3 = RS-232 interface

The sensor type response is a 0 for a request to an M15

4.2.12. ASCII m Request for OEM version number (complete)

mS

This command is used for Medicaid testing, not for general use.

4.2.13. ASCII n Request for self test (incomplete, x = data, chk = checksum)

n x chk

| | | _____ modulo 32 checksum

| | _____ test type

| _____ "self test" character

The test types available from 0 to 31.

0 = Lo Batt cal

1 = Display test

2 = EEPROM test

3 = Initialization of EEPROM

Revision History

Rev 0.6 Added baud rate modifier

Rev 0.7 Added command "m" for requesting OEM connected rev

Added command "n" for requesting self testes, cannot be used by OEM

Rev 0.8 Added System Failure and No module Attached to error codes

Rev 0.9 Added error code for analog output calibration failure, used on M15As

Added note in sign-off sheet for electronic copy

Rev 1.0 Removed Temperature and Spare output, changed the checksum offset.

Added the perfusion send interval time from 1 to 255.

Removed request of serial number.

Incorporated Perfusion Quality values range, 1-10.

Table 1: Perfusion Quality values

Pulse Amplitude	Perfusion Quality
0.0% -0.9%	0-6
1.0%-2.0%	7
3.0%-4.0%	8
5.0%-10.0%	9
10%-20%	10