

## THM010-HT001

Thermometer, Non-Contact Infrared Human



## Key Features:

- Precision: Measuring tolerance is less than  $\pm 0.3$  °C
- Swiftness: Measuring times is less than 1 second.
- Easy to operate: One key measurement, operate conveniently
- Non-contact: Forehead measurement, no skin contact.
- Large screen display : Large LCD screen display. Makes clear digital performance in any ray.
- High temperature alarm
- Data storage: save 32 group of measuring data, easy for analysis and comparison
- Three colour back lights for different temperature zone.
- Language available : English, French, Spanish.

## **Product Specifications**

Normal conditions of use

Circumstance temperature:10°C ~ 40°C (50°F ~ 104°F) Humidity range:  $\leq$  85%

Power: ≤300mW

**Size:** 175 x 130 x 65mm (L x W x H)

Weight: 110g

Display Resolution: 0.1°C (0.14°F)

Measuring range:

In body mode: 32°C - 42.9°C (89.6°F-109.2°F)
In Surface Temp mode: 0°C - 60°C (32°F - 140°F)

Precision: From 32°C to 36°C ±0.3°C 36°C to 39°C ±0.2°C

36°C to 39°C ±0.2°C 39°C to 42.9°C ±0.3°C

Consumption: < 300mW

Measuring distance: 3-5cm (1.2in - 2in)

Automatic stop: 5 sec



17517 Fabrica Way, Suite H, Cerritos, CA 90703, USA (Tel) 714-367-2848 (Fax) 714-367-2852 www.mediaidinc.com info@mediaidinc.com