

# MODEL 100 SERIES

## Mini Pulse Oximeters



Model 100  
with Integral  
Finger Sensor

### Two Configurations - Cordless or Cabled -----

Series 100 oximeters can be used with either the Integral Sensors or with the Cable Adapter Module and a wide variety of adult and pediatric sensors. A patented feature permits an easy interchange between the two sensor configurations. Opening the Integral Sensor for finger placement automatically turns on the device. The large red LED display alternates between oxygen saturation and pulse rate every few seconds.

### Model 100 -----

The Model 100 is the basic spot checking mini oximeter of the series and can be relied upon for fast and accurate measurement. A single "AA" alkaline battery gives approximately 1200 spot checks (Duracell Ultra recommended).

### Model 130 -----

The Model 130 has memory to retain oxygen saturation and pulse rate data, and can operate in three different modes to best suit your needs:

1. The Model 130 can retain up to 30 minutes of information for a single patient, making the device ideally suited for treadmill or other long tests.
2. The Model 130 can also be used in walk around tests, sorting data for up to 20 consecutive patients, tagging each reading with a patient number, time and date.
3. The Model 130 also includes expanded memory for sleep screening where data can be stored every 5 seconds for upto 10.5 hours.

Data stored in any mode can be printed or quickly transferred by infrared to a personal computer.



Model 130  
with cable  
Adapter and  
Soft Sensor

## Equipment Specifications

### Oxygen Saturation & Pulse

Display	Three digit, seven segment LED Alternates between Oxygen Saturation and Pulse Rate
Resolution	SpO <sub>2</sub> 1% Pulse 1 BMP
Range	SpO <sub>2</sub> 20-100% Pulse 25-250 BMP
Accuracy	SpO <sub>2</sub> 100 - 70%, ± 2% 69 - 60%, ±3% less than 60% unspecified Pulse 25 - 200 BPM, ±2BPM or 2% (which ever is greater) Above 200, ±3%
Response Time	SpO <sub>2</sub> 8 sec for 80% of patients Pulse 8 sec for 80% of patients
Update Time	Continuous

### Demensions

Size	4.7(L) x 1.8(W) x 1.0(H) Inch 12(L) x 4.7(W) x 2.5(H) cm
Weight	3.45 oz, 98 gm with either Integral Finger Sensor or Cable Adapter Module (weight includes battery)

### Power Source

Battery Type	one 1.5V "AA" Alkaline battery
Battery Life	Approximately 1200 spot checks

### Warranty

Model 100 & 130 each carry a two year warranty.  
(sensor and accessories not included)



tabular format with relevant statistics. A zoom feature allows analysis in detail and comments may be added to the permanent file.

### OxySoft™ Software -----

OxySoft™ pulse oximetry software is available on a CD-ROM and has been developed by Mediaid Inc. to display saturation and pulse rate data in a variety of helpful formats. OxySoft™ software enables the health care professional to present data in graphical or

17517 Fabrica Way Suite H Cerritos, CA 90703 USA  
(Tel) 714-367-2848 (Fax) 714-367-2852  
[www.mediaidinc.com](http://www.mediaidinc.com) [info@mediaidinc.com](mailto:info@mediaidinc.com)

