

MEDIAID

# THM010-VT001

Thermometer, Non-Contact  
Infrared Veterinary



## Key Features:

- Precision: Measuring tolerance is less than  $\pm 0.3^{\circ}\text{C}$
- Swiftness: Measuring times is less than 1 second.
- Easy to operate: One key measurement, operate conveniently
- Non-contact: No skin contact.
- Large screen display : Large LCD screen display.  
white back-light, Makes clear digital performance in any ray.
- High temperature alarm
- Language available : English, French, Spanish.
- Memory of last 32 readings.



## Product Specifications

### Normal conditions of use

Operating temperature:  $10^{\circ}\text{C} \sim 40^{\circ}\text{C}$  ( $50^{\circ}\text{F} \sim 104^{\circ}\text{F}$ )  
Humidity range:  $\leq 85\%$

**Power:** DC 3V (2batteries AAA)

**Size:** 141 x 44 x 52mm - 5.6 x 1.4 x 2.1in (L x W x H)

**Weight:** 110g

**Display Resolution:**  $0.1^{\circ}\text{C}$  ( $0.14^{\circ}\text{F}$ )

### Measuring range:

In body mode :  $32^{\circ}\text{C} - 42.9^{\circ}\text{C}$  ( $89.6^{\circ}\text{F} - 109.2^{\circ}\text{F}$ )  
In Surface Temp mode :  $0^{\circ}\text{C} - 60^{\circ}\text{C}$  ( $32^{\circ}\text{F} - 140^{\circ}\text{F}$ )

**Precision :** From  $36^{\circ}\text{C}$  to  $39^{\circ}\text{C}$  ( $96.8^{\circ}\text{F} - 102.2^{\circ}\text{F}$ ) =  $\pm 1.0^{\circ}\text{C}$  ( $0.4^{\circ}\text{F}$ )

**Consumption :**  $\leq 300\text{mW}$

**Accuracy :**  $\pm 0.3^{\circ}\text{C}$  ( $0.6^{\circ}\text{F}$ )

**Measuring distance :** 3-5cm (1.2in - 2in)

**Automatic stop :** 5 sec

### Mediaid Inc.

17517 Fabrica Way, Suite H,  
Cerritos, CA 90703, USA

(Tel) 714-367-2848 (Fax) 714-367-2852

[www.mediaidinc.com](http://www.mediaidinc.com) [info@mediaidinc.com](mailto:info@mediaidinc.com)