

MEDIAID

Veterinary Pulse Oximeter Sensors

Lingual Clip Sensor

Mediaid's Lingual Clip Sensor is a reusable veterinary sensor suitable for all dogs, cats, and various large animals. The Replaceable silicon pads on this sensor offer added grip to the tongue, which limits sensor movement and reduces interruptions during monitoring. The sensor is packaged with six extra silicon pads. The lingual clip sensor and silicon pads are water resistant and can be easily cleaned. This sensor provides only oximetry data.

- RJ12 Connector
REF/PN POX022-450
(use with models 340V, 5340V)
REF/PN POX022-455
(use with models 500V, 500VP)

Silicon Pads

The Mediaid Silicon Pads for the lingual clip are removable pads that can be used multiple times. These pads provide additional adherence to the tongue to limit interruptions while monitoring. The package contains 24 water resistant pads.

- REF/PN POX025-450

Rectal Sensor

This veterinary rectal sensor provides pulse oximetry and core temperature data. It is suitable for both large and small animals. A thermistor in the sensor allows body temperature to be displayed on the pulse oximeter. Temperature can be displayed intermittently every two minutes or on demand (°C or °F). A double beep sounds each time the temperature is displayed.

- RJ12 Connector
REF/PN POX022-650
(use with models 340V, 5340V)
REF/PN POX022-655
(use with models 500V, 500VP)

Dual Sensor Adapter

The Dual Sensor Adapter is compatible with a Mediaid veterinary pulse oximeter. This dual sensor enables the rectal sensor and lingual clip sensor to be used simultaneously, providing core temperature and pulse oximetry data for animals where feces may obstruct oximetry readings in the rectal probe.

- RJ12 Connector
REF/PN POX025-700

17517 Fabrica Way, Suite H; Cerritos, CA 90703 USA

p 714.367.2848 f 714.367.2852

www.mediaidinc.com info@mediaidinc.com