



# Veterinary Compatible Sensors Chart













Compatibility	Monitors	OEM Part Number	Lingual Sensor	Rectal Sensor	Cable Length in Feet	Connector Type		
								
<b>NELLCOR</b>	N-20PAV, N-65VSAT, N-560VSAT, NPV-70, N-85	VSAT (Lingual Clip Sensor)	CST090-51209	CST080-51209	9 feet (3 meters)	DB9		
<b>NONIN</b>	9847V, 8500AV, LifeSenseVet, RespSenseVet, PulseSenseVet	2000SL (Lingual Clip Sensor)	CST090-51259	CST080-51259	9 feet (3 meters)	DB9		
<b>SURGIVET</b>	V3395, V3402, V3304, V3404, V3204P, V60041, V60046, V90041, V90043, V9200, V9201, V9203, V9204	V1700 (Reflectance Sensor)V1703 (Universal "Y" Lingual Sensor)	CST090-51049	CST080-51049	9 feet (3 meters)	DB9		
<b>MINDRAY with Mindray Technology</b>	MEC-1200Vet, PM-9000Vet, PM-60Vet, PM-50Vet	0010-10-43237	CST090-51409	CST080-51409	9 feet (3 meters)	REDEL 6 PIN 40 DEG DOUBLE NOTCH		
<b>MINDRAY with Nellcor Technology</b>	MEC-1000, PM-9000	Not available	CST090-51209	CST080-51209	9 feet (3 meters)	DB9		
<b>GOLDWAY with BCI Technology</b>	VET-400A, VET-600A, VET 420F, VET 400C	3078 (BCI Small Lingual Sensor) V-3222 (BCI Large Lingual Sensor)	CST090-51309	CST080-51309	9 feet (3 meters)	DB9		

MEDIAID INC., 17517 Fabrica Way Suite H, Cerritos, CA 90703, USA. Tel: 714-367-2848, Fax: 714-367-2852, Email: info@mediaidinc.com, url: www.mediaidinc.com

# Veterinary Compatible Sensors Chart



Compatibility	Monitors	OEM Part Number	Lingual Sensor	Rectal Sensor	Cable Length in Feet	Connector Type		
								
<b>GOLDWAY with Nellcor Technology</b>	VET-400A, VET-600A, VET 420F, VET 400C	VSAT (Lingual Clip Sensor)	CST090-51209	CST080-51209	9 feet (3 meters)	DB9		
<b>CRITICARE</b>	503 Spot, 504 Plus, 504DX, 506DXN, 501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), 8100, Scholar	Not available	CST090-51079	CST080-51079	9 feet (3 meters)	DB9		
<b>HESKA</b>	4400, 4402, 4404, 4700	CAT4060, CAT0451	CST090-51609	CST080-51609	9 feet (3 meters)	DB9		
<b>BIONET</b>	BM5 VET, BM3 VetNext, BM3 Vet, BM3 VetPlus, OXY9 Vet	B-SPO2-V	CST090-51809	CST080-51809	9 feet (3 meters)	DB9	