






International Cross Reference Chart for Reusable Sensors with Mediad Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
ARMSTRONG	AD-900, 1000, 2000, 3000, 4000	BCI 3044	 CST060-3104	 CST062-3104	 CST064-3104	 CST070-3104	 CST100-3104	 CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
ARTEMA	MM200		 CST060-31019	 CST062-31019	 CST064-31019	 CST070-31019	 CST100-31019	 CST102-31019	9 / 270	NICOLAY 12 PIN		
ARTEMA	Diascope	Nelcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
BAXTER	ASAT	Simed 600015	 CST060-31029	 CST062-31029	 CST064-31029	 CST070-31029	 CST100-31029	 CST102-31029	9 / 270	HYPERTRONICS 9 PIN		
BAXTER	Explorer, Monarch	Nelcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
BCI	3100, 3101, 6100, 9100, Advisor, Autocorr, Minitorr, Minitorr Plus		 CST060-31039	 CST062-31039	 CST064-31039	 CST070-31039	 CST100-31039	 CST102-31039	9 / 270	REDEL 7 PIN SINGLE NOTCH		
BCI	3100, 3101, 3300, 3301, 3302, 3303, 3304, 3401, 6100, 9100, Advisor, Autocorr, Minitorr, Minitorr Plus	BCI 3044	 CST060-3104	 CST062-3104	 CST064-3104	 CST070-3104	 CST100-3104	 CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
BCI	Monitors using BCI Module		 CST060-31059	 CST062-31059	 CST064-31059	 CST070-31059	 CST100-31059	 CST102-31059	9 / 270	REDEL 7 PIN 40 DEG DOUBLE NOTCH		
BCI	Monitors using BCI Module		 CST060-31089	 CST062-31089	 CST064-31089	 CST070-31089	 CST100-31089	 CST102-31089	9 / 270	REDEL 5 PIN SINGLE NOTCH		





















International Cross Reference Chart for Reusable Sensors with Mediad Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
BCI	Monitors using BCI Module		 CST060-31559	 CST062-31559	 CST064-31559	 CST070-31559	 CST100-31559	 CST102-31559	9 / 270	REDEL 7 PIN 80 DEG DOUBLE NOTCH		
BIOSYS	BPM-200	Nellcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			 CST060-31209	 CST062-31209	 CST064-31209	 CST070-31209	 CST100-31209	 CST102-31209	9 / 270			
BIRD	4400	CSI 511-10L	 CST060-31069	 CST062-31069	 CST064-31069	 CST070-31069	 CST100-31069	 CST102-31069	9 / 270	LEMO 5 PIN		
BIRD	4400	CSI 511SD	 CST060-3107	 CST062-3107	 CST064-3107	 CST070-3107	 CST100-3107	 CST102-3107	3 / 90	DB9		
			 CST060-31079	 CST062-31079	 CST064-31079	 CST070-31079	 CST100-31079	 CST102-31079	9 / 270			
BPL	OM 5509		 CST060-31059	 CST062-31059	 CST064-31059	 CST070-31059	 CST100-31059	 CST102-31059	9 / 270	REDEL 7 PIN 40 DEG DOUBLE NOTCH		
BPL	MPM 5633		 CST060-31089	 CST062-31089	 CST064-31089	 CST070-31089	 CST100-31089	 CST102-31089	9 / 270	REDEL 5 PIN SINGLE NOTCH		
BPL	OM 5609, PO 5529	Ohmeda 6051-0000-035	 CST060-31309	 CST062-31309	 CST064-31309	 CST070-31309	 CST100-31309	 CST102-31309	9 / 270	HYPERTRONICS 9 PIN		
CAS MEDICAL	OscilloMate 9002, 9303, 9302S	Nonin 8000AA-1	 CST060-3125	 CST062-3125	 CST064-3125	 CST070-3125	 CST100-3125	 CST102-3125	3 / 90	DB9		
			 CST060-31259	 CST062-31259	 CST064-31259	 CST070-31259	 CST100-31259	 CST102-31259	9 / 270			
COLIN	Press-Mate SAT, BP 88, 306, 508, 808, 88S	Nellcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			 CST060-31209	 CST062-31209	 CST064-31209	 CST070-31209	 CST100-31209	 CST102-31209	9 / 270			




International Cross Reference Chart for Reusable Sensors with Mediaid Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
COROMETRICS	Corometrics 118, 120, 128, 128F, 129, 129F, 188, 510, 511	Nellcor DS100A							3 / 90	DB9		
			CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	9 / 270			
CRITICARE	501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), 8100, Scholar	CSI 511-10L							9 / 270	LEMO 5 PIN		
CRITICARE	503 Spot, 504Plus, 504DX, 506DXN, 501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), 8100, Scholar	CSI 511SD							3 / 90	DB9		
			CST060-3107	CST062-3107	CST064-3107	CST070-3107	CST100-3107	CST102-3107	9 / 270			
CRITIKON	Critikon Dinamap 9700, 9710, 9720 Plus, Compact, MPS, Select, Dinamap Pro 100-400	Nellcor DS100A							3 / 90	DB9		
			CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	9 / 270			
CRITIKON	Critikon Dinamap 8700, 8710, 8720 Plus, Oxyshuttle	Critikon 9084							3 / 90	COMPUSHIELD		
DATASCOPE	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascop SpO ₂), Passport (Datascop SpO ₂)	Datascop 0600-00-0026-01							3 / 90	DB9		
			CST060-3109	CST062-3109	CST064-3109	CST070-3109	CST100-3109	CST102-3109	9 / 270			
DATASCOPE	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascop SpO ₂), Passport (Datascop SpO ₂)	Datascop 0020-00-071-01							9 / 270	DIN 8		
DATASCOPE	Expert, Passport (Nellcor SpO ₂), Accutorr Plus (Nellcor SpO ₂)	Nellcor DS100A							3 / 90	DB9		
			CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	9 / 270			
DATEX	AS/3, CS/3, CH/S, CH/RS, CG/2S, Cardiacap II, Capnomac, Ultima, Light Monitor, M-Estr, Oscar II, Satellite Trans Plus, ULTSVI, ULTSVO, ULTSIO, ULTS00, OSP-200	Datex SAS-F4							9 / 270	DPK1		

International Cross Reference Chart for Reusable Sensors with Mediaid Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	DB9	Connector Type	
												
DATEX	AS/3, CS/3, CH/S, CH/RS, CG/2S, Cardicap II, Capnomac, Ultima, Light Monitor, M-Estpr, Oscar II, Satellite Trans Plus, ULTSVI, ULTSVO, ULTSIO, ULTS00,OSP-200	Datex SAS-F	CST060-3114	CST062-3114	CST064-3114	CST070-3114	CST100-3114	CST102-3114	3 / 90	DB9		
			CST060-31149	CST062-31149	CST064-31149	CST070-31149	CST100-31149	CST102-31149	9 / 270			
DIGICARE	Digimax 5500	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
DIMA ITALIA	All Models	Nelcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
DRAGER	Drager Oxipac 2000	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
DRAGER	Vitalert 1000, Cicero, Narkomed	Nelcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
ELMED	PSO 1000	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
ENTRI MEDICAL	Propaq 102, 104, 106	Nelcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
FUKUDA DENSHI	HG302 Module, DS-5100E, DS-5300W	Nelcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
GAMBRO	BP 806	Nelcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			


International Cross Reference Chart for Reusable Sensors with Mediaid Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
												
GOLDWAY	UT4000A	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
HEALTHDYNE	All Models	Nonin 8000AA-1	CST060-3125	CST062-3125	CST064-3125	CST070-3125	CST100-3125	CST102-3125	3 / 90	DB9		
			CST060-31259	CST062-31259	CST064-31259	CST070-31259	CST100-31259	CST102-31259	9 / 270			
INNOMED	Oxycard, Innocare-T0	Nonin 8000AA-1	CST060-3125	CST062-3125	CST064-3125	CST070-3125	CST100-3125	CST102-3125	3 / 90	DB9		
			CST060-31259	CST062-31259	CST064-31259	CST070-31259	CST100-31259	CST102-31259	9 / 270			
INVACARE	All Models	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
INVIVO RESEARCH	4500 Plus 2	Nelcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
INVIVO RESEARCH	4500, 4500 Plus, Millenia, Sat Scout Omega 1445, Omni-Trak TVS	Invivo 9383	CST060-31419	CST062-31419	CST064-31419	CST070-31419	CST100-31419	CST102-31419	9 / 270	LEMO 7 PIN		
IVY BIOMEDICAL	405A, 405D, 405P, 405T	Nelcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
L & T	All Models		CST060-31219	CST062-31219	CST064-31219	CST070-31219	CST100-31219	CST102-31219	9 / 270	REDEL 7 PIN 80 DEG DOUBLE NOTCH		
LYNN MEDICAL	All Models	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			

International Cross Reference Chart for Reusable Sensors with Mediad Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
MAGNA MEDICAL	MON-MTEK-30, 33, 35, 40, 85, 1200, 1201, 1300, 1401, 1406, 1411, 1413, 1501, 1503, MON-MTEK-MOD2, MOD3	BCI 3044	 CST060-3104	 CST062-3104	 CST064-3104	 CST070-3104	 CST100-3104	 CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
MARQUETTE EAGLE	Eagle, Dash, SOLAR, TRAM Modular 250SL, 450SL, 650SL, 850A, 850SL	Nelcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
MARQUETTE EAGLE	Eagle, Dash, SOLAR, TRAM Modular 250SL, 450SL, 650SL, 850A, 850SL		 CST060-31359	 CST062-31359	 CST064-31359	 CST070-31359	 CST100-31359	 CST102-31359	9 / 270	EAGLE FLAT 11 PIN		
MARQUETTE HELLIGE	Hellige Eagle 1000, .3000, SMU, SMK, SMV	Nelcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
MARQUETTE HELLIGE	7000 Series, TRAM A.R., TRAM 200SL, 400SL, 600SL, 800A, 800SL		 CST060-31329	 CST062-31329	 CST064-31329	 CST070-31329	 CST100-31329	 CST102-31329	9 / 270	HYPERTRONICS 9 PIN		
MARQUETTE HELLIGE	Hellige Eagle 1000, 3000, SMU, SMK, SMV		 CST060-31459	 CST062-31459	 CST064-31459	 CST070-31459	 CST100-31459	 CST102-31459	9 / 270	HELLIGE ROUND 10 PIN		
M.D.E.	Escort 100 OPT11A, Escort 300A, Escort II OPT30, Prism	Nelcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
MEDIANA	YM-2002 Series	Nelcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
MEDLAB	NANOX2	BCI 3044	 CST060-3104	 CST062-3104	 CST064-3104	 CST070-3104	 CST100-3104	 CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			





















International Cross Reference Chart for Reusable Sensors with Medicaid Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
MEDLAB	POX 10								9 / 270	REDEL 7 PIN 80 DEG DOUBLE NOTCH		
MEK	MP-100, MP-110, MP-400, MP-1000		CST060-31559	CST062-31559	CST064-31559	CST070-31559	CST100-31559	CST102-31559	3 / 90	DB9		
			CST060-31519	CST062-31519	CST064-31519	CST070-31519	CST100-31519	CST102-31519	9 / 270			
MEK	MP-100 Plus, MP-111		CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
MENNEN	Horizon XL (Nellcor SpO ₂), Mercury 1200, MR 1303	Nellcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
MENNEN	Horizon XL (Ohmeda SpO ₂)	Mennen 0800-060-540	CST060-31309	CST062-31829	CST064-31829	CST070-31829	CST100-31829	CST102-31829	9 / 270	HYPERTRONICS 9 PIN		
MINDRAY	MEC-1000, PM-9000		CST060-3162	CST062-31	CST064-31	CST070-31	CST100-31	CST102-31	9 / 270	REDEL 6 PIN 40 DEG DOUBLE NOTCH		
NELLCOR	N20, 100, 180, 190, 200, 390, 3000, 6000, NPB-290, 295	Nellcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
NELLCOR	Monitors using Nellcor Module		CST060-31219	CST062-31219	CST064-31219	CST070-31219	CST100-31219	CST102-31219	9 / 270	REDEL 7 PIN 80 DEG DOUBLE NOTCH		
NIHON KOHDEN	AL-900P Module, BSM-1102, 9510, OGS-2000, BSM-4104A, 4114A	Nellcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			



























International Cross Reference Chart for Reusable Sensors with Mediaid Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
NIHON KOHDEN	AL-800PA Module, BSM-1101, 2101/2, 3101, 7100, 7200, OLV-1100A, Lifescope	Nihon Kohden TL-101S							9 / 270	DPK1 10 PIN		
NIHON KOHDEN	AL-800PA Module, BSM-1101, 2101/2, 3101, 7100, 7200, OLV-1100A, Lifescope		CST060-3124	CST062-3124	CST064-3124	CST070-3124	CST100-3124	CST102-3124	3 / 90	DB9		
			CST060-31249	CST062-31249	CST064-31249	CST070-31249	CST100-31249	CST102-31249	9 / 270			
NONIN	8500, 8600, 8700, 8800	Nonin 8000AA-1	CST060-3125	CST062-3125	CST064-3125	CST070-3125	CST100-3125	CST102-3125	3 / 90	DB9		
			CST060-31259	CST062-31259	CST064-31259	CST070-31259	CST100-31259	CST102-31259	9 / 270			
NOVAMETRIX	505, 510, 511, 515 B/C, 520, 7000, 7100, Oxypleth	Novamatrix 8776-00	CST060-31269	CST062-31269	CST064-31269	CST070-31269	CST100-31269	CST102-31269	9 / 270	HYPERTRONICS 7 PIN		
NOVAMETRIX	505, 510, 511, 515 B/C, 520, 7000, 7100, Oxypleth		CST060-3128	CST062-3128	CST064-3128	CST070-3128	CST100-3128	CST102-3128	3 / 90	DB9		
			CST060-31289	CST062-31289	CST064-31289	CST070-31289	CST100-31289	CST102-31289	9 / 270			
ODAM	All Models	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
OHMEDA	3700, 3710, 3770, 3775, Biox 3740, 3760, 3800, 3900, 4700, Tuffsat, OxiCap, 5250 RGM, Rascal II	Ohmeda 6051-0000-035	CST060-3130	CST062-3130	CST064-31309	CST070-3130	CST100-3130	CST102-3130	3 / 90	HYPERTRONICS 9 PIN		
			CST060-31309	CST062-31309		CST070-31309	CST100-31309	CST102-31309	9 / 270			
OHMEDA	3700, 3710, 3770, 3775, Biox 3740, 3760, 3800, 3900, 4700, Tuffsat, OxiCap, 5250 RGM, Rascal II		CST060-3131	CST062-3131	CST064-3131	CST070-3131	CST100-3131	CST102-3131	3 / 90	DB9		
			CST060-31319	CST062-31319	CST064-31319	CST070-31319	CST100-31319	CST102-31319	9 / 270			
PACETECH	Vitalmax 800 Series (Nellcor SpO ₂)	Nellcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			

International Cross Reference Chart for Reusable Sensors with Medicaid Part Numbers



Make	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Ear Sensor	Y - Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet / Cm.	Connector Type		
PHYSIO-CONTROL	Lifepak 12	Nellcor DS100A	 CST060-3120	 CST062-3120	 CST064-3120	 CST070-3120	 CST100-3120	 CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
RGB MEDICAL	Omicrom FT, Omicrom PREMIUM		CST060-31039	CST062-31039	CST064-31039	CST070-31039	CST100-31039	CST102-31039	9 / 270	REDEL 7 PIN SINGLE NOTCH		
RGB MEDICAL	Omicrom FT, Omicrom PREMIUM	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
SCHILLER	CM8, Argus TM-7	Nellcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
SIEMENS	Sirec 700, 630, Sirem SpO ₂ Module, SC 6002XL, 7000, 8000, 9000XL	Nellcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			
SIMED	S-50, S-100, S-100e	Simed 600015	CST060-31029	CST062-31029	CST064-31029	CST070-31029	CST100-31029	CST102-31029	9 / 270	HYPERTRONICS 9 PIN		
SPACELABS BURDICK	PRO ₂	BCI 3044	CST060-3104	CST062-3104	CST064-3104	CST070-3104	CST100-3104	CST102-3104	3 / 90	DB9		
			CST060-31049	CST062-31049	CST064-31049	CST070-31049	CST100-31049	CST102-31049	9 / 270			
SPACELABS	90419-01, 90651A-07, 90351-D/G, 90465, 90466, 90467, 90489, 90651-07/08, 90651A, 90496 Ultraview (Novamatrix SpO ₂)	Spacelabs 015-0130-00	CST060-31269	CST062-31269	CST064-31269	CST070-31269	CST100-31269	CST102-31269	9 / 270	HYPERTRONICS 7 PIN		
SPACELABS	90419-01, 90651A-07		CST060-3128	CST062-3128	CST064-3128	CST070-3128	CST100-3128	CST102-3128	3 / 90	DB9		
			CST060-31289	CST062-31289	CST064-31289	CST070-31289	CST100-31289	CST102-31289	9 / 270			
WELCH ALLYN	Propaq 206-EL, 204-EL, 202-EL, 106-EL	Nellcor DS100A	CST060-3120	CST062-3120	CST064-3120	CST070-3120	CST100-3120	CST102-3120	3 / 90	DB9		
			CST060-31209	CST062-31209	CST064-31209	CST070-31209	CST100-31209	CST102-31209	9 / 270			