
















USA CROSS REFERENCE CHART FOR REUSABLE SENSORS WITH MEDIAID PART NUMBERS



Compatibility	Monitors	Adult Finger Sensor OEM Part Number	Adult Finger Sensor	Pedi Finger Sensor	Y Sensor	Soft Adult Sensor	Soft Pedi Sensor	Cable Length in Feet	Connector Type		
											
BCI	3100, 3101, 3300, 3301, 3302, 3303, 3304, 3401, 6100, 9100, Advisor, Autocorr, Minitorr, Minitorr Plus	BCI 3044	CST060-3104	CST062-3104	CST070-3104	CST100-3104	CST102-3104	3	DB9		
			CST060-31049	CST062-31049	CST070-31049	CST100-31049	CST102-31049	9			
NELLCOR®	N20, 100, 180, 190, 200, 390, 3000, 6000, NBP-290, 295	NELLCOR DS100A	CST060-3120	CST062-3120	CST070-3120	CST100-3120	CST102-3120	3	DB9		
			CST060-31209	CST062-31209	CST070-31209	CST100-31209	CST102-31209	9			
NONIN	8500, 8600, 8700, 8800	Nonin 8000AA-1	CST060-3125	CST062-3125	CST070-3125	CST100-3125	CST102-3125	3	DB9		
			CST060-31259	CST062-31259	CST070-31259	CST100-31259	CST102-31259	9			
OHMEDA	3700, 3710, 3700, 3775, Blox 3740, 3760, 3800, 3900, 4700, Tuffsat, OxiCap, 5250 RGM, Rascal II	Ohmeda 6054-0000-035	CST060-3130	CST062-3130	CST070-3130	CST100-3130	CST102-3130	3	HYPERTRONICS 9 PIN		
			CST060-31309	CST062-31309	CST070-31309	CST100-31309	CST102-31309	9			
OHMEDA	3700, 3710, 3700, 3775, Blox 3740, 3760, 3800, 3900, 4700, Tuffsat, OxiCap, 5250 RGM, Rascal II		CST060-3131	CST062-3131	CST070-3131	CST100-3131	CST102-3131	3	DB9		
			CST060-31319	CST062-31319	CST070-31319	CST100-31319	CST102-31319	9			

MEDIAID INC., 17517 Fabrica Way Suite H, Cerritos, CA 90703, USA. Tel: 714-367-2848, Fax: 714-367-2852, Email: info@mediaidinc.com, url: www.mediaidinc.com