

Infrared Forehead Thermometer Non-Contact Measurement

THM010-HT002

- Illuminated Display
- Quick & Accurate
- Easy to Use
- 10 Memories



FEATURES

Display: illuminated display

Fast: measuring time ≤ 1 s

Easy operation: simple one button press

Energy saving: automatically power off in 1 minute

Long life: use over 2500 times with 1PC of CR2032 button battery

Dual temperature mode: change between C and F

Memory function: the last 10 memories

Beeper signal: on/off, test result and low voltage warning

Visual warning: low voltage warning

Self-test: power on, all segments will display on the LCD in 1s.

Non-contact measurement: no touch forehead measuring

Can be used for: all ages

Box content: 1 pc Infrared Forehead thermometer, instruction sheet, storage box

SPECIFICATIONS

Operating ambient:
Temperature: $10^{\circ}\text{C} \sim 40^{\circ}\text{C}$
($50.0^{\circ}\text{F} \sim 104.0^{\circ}\text{F}$)
Humidity: $\leq 85\% \text{RH}$

The display range:
 $32.0^{\circ}\text{C} \sim 42.9^{\circ}\text{C}$
($89.6^{\circ}\text{F} \sim 109.2^{\circ}\text{F}$)

Accuracy:
 $\pm 0.3^{\circ}\text{C}$ ($32.0^{\circ}\text{C} \sim 42.9^{\circ}\text{C}$)
 $\pm 0.5^{\circ}\text{F}$ ($89.6^{\circ}\text{F} \sim 109.2^{\circ}\text{F}$)

Storage and transport condition:
temperature: $-20^{\circ}\text{C} \sim 55^{\circ}\text{C}$
($-4^{\circ}\text{F} \sim 131^{\circ}\text{F}$)
Humidity: $\leq 93\% \text{RH}$

Display resolution: 0.1°C (0.1°F)

Power consumption: $\leq 30\text{mW}$

Battery: 1PCCR2032 button battery

Dimension: $4.13 \times 1.38 \times 1.38$ inch
($105 \times 3 \times 35$ mm)

Net weight: 1.3 oz (36g)

Temperature unit: $^{\circ}\text{C}$ or $^{\circ}\text{F}$



Faichney Medical
17517 Fabrica Way
Cerritos, CA 90703

www.Faichney-Medical.com
Info@Faichney-Medical.com